

# Delay Alternating with Nutation for Tailored Excitation (DANTE)

The DANTE-prepared multicontrast 2D [black blood imaging](#) technique is a promising new tool for [MR](#) imaging of [carotid artery](#) walls. Additionally, the 3D DASH sequence enables 3D high-spatial-resolution fast T1-weighted imaging of carotid artery walls.

From:

<https://neurosurgerywiki.com/wiki/> - **Neurosurgery Wiki**

Permanent link:

[https://neurosurgerywiki.com/wiki/doku.php?id=alternating\\_with\\_nutation\\_for\\_tailored\\_excitation](https://neurosurgerywiki.com/wiki/doku.php?id=alternating_with_nutation_for_tailored_excitation)

Last update: **2024/06/07 02:54**

