

CD HORIZON LEGACY SPINAL SYSTEM

The CD Horizon® Legacy™ Spinal System can be used to treat [scoliosis](#).



SIZE AND PLACEMENT

The CD Horizon Legacy Spinal System consists of [rods](#), [hooks](#), and [screws](#) for implantation in the spine to correct the abnormal curvature. The system is available in multiple rod diameters and screw sizes so your surgeon can choose the appropriate size for your condition, anatomy, and activity level.

The CD Horizon Legacy Spinal System has been used to treat scoliosis since 2004. Till 2016 there have been over 500,000 surgeries performed with the CD Horizon Legacy Spinal System.

The CD Horizon Legacy Spinal System is made out of [titanium](#) or stainless steel implantable grade metal.

From:

<https://neurosurgerywiki.com/wiki/> - Neurosurgery Wiki

Permanent link:

https://neurosurgerywiki.com/wiki/doku.php?id=cd_horizon

Last update: **2024/06/07 02:53**

