

COSEAL

<http://www.coseal.com/>

Translucent seal Helps to keep suture lines visible.

Can be sutured through

Can be sutured through After polymerization, additional sutures can be placed through it without disturbing the seal.¹

No refrigeration required

No refrigeration required Can be stored at room temperature (approximately 25 °C).²

Multiple application options

Multiple application options May be used for focal and broad applications.² Use caution when applying with pressurized gas.²

From:

<https://neurosurgerywiki.com/wiki/> - **Neurosurgery Wiki**

Permanent link:

<https://neurosurgerywiki.com/wiki/doku.php?id=coseal>

Last update: **2024/09/18 06:44**

