

The Dual-Luciferase® Reporter (DLR™) Assay System provides an efficient means of performing two reporter assays. ... After quantifying the firefly luminescence, this reaction is quenched, and the Renilla luciferase reaction is initiated simultaneously by adding Stop & Glo® Reagent to the same sample.

From:

<https://neurosurgerywiki.com/wiki/> - **Neurosurgery Wiki**

Permanent link:

https://neurosurgerywiki.com/wiki/doku.php?id=dual-luciferase_reporter_assay

Last update: **2024/06/07 02:49**

