

ETurn® TLIF

<https://www.icotec-medical.com>

The BlackArmor® material is radiolucent and rapidly integrates with the vertebral endplates thanks to the porous Ti-iT® titanium coating.

The icotec ETurn® TLIF Cage may be implanted via the surgeon's preferred posterior approach: [TLIF](#) (transforaminal lumbar interbody fusion), [EPLIF](#) (extraforaminal posterior lumbar interbody fusion) or [PLIF](#) (posterior lumbar interbody fusion) using biportal cage implantations.

From:

<https://neurosurgerywiki.com/wiki/> - **Neurosurgery Wiki**

Permanent link:

<https://neurosurgerywiki.com/wiki/doku.php?id=eturn>

Last update: **2024/06/07 02:56**

