

Markham-Meyerding Hemilaminectomy Retractor



This retractor is also commonly used in lumbar discectomy. It consists of two **blades**, a solid blade that is inserted beneath the **transverse process** and a hollow blade that is inserted beneath the superior facet.

From:
<https://neurosurgerywiki.com/wiki/> - Neurosurgery Wiki

Permanent link:
https://neurosurgerywiki.com/wiki/doku.php?id=markham-meyerding_hemilaminectomy_retractor

Last update: 2024/06/07 02:57

