

Pedicle Screw Preparation

[Pedicle Screw](#) Preparation is performed utilizing a selection of [Awls](#), [Pedicle Probes](#), Ball Tip Feelers and [Bone Taps](#).

- Probes and Bone Taps are marked to indicate the appropriate length of [Polyaxial Screw](#).

From:

<https://neurosurgerywiki.com/wiki/> - **Neurosurgery Wiki**

Permanent link:

https://neurosurgerywiki.com/wiki/doku.php?id=pedicle_screw_preparation

Last update: **2024/06/07 02:55**

