

# Pulsed arterial spin labeling

In pulse [arterial spin labeling](#) (PASL), blood water is inverted as it passes through a labeling slab (of 15 to 20 cm) instead of a plane. There are different variations of this implementations, including EPISTAR and PICORE and PULSAR. Most scanner have been designed to have PASL out-of-the-box for research use.

From:

<https://neurosurgerywiki.com/wiki/> - **Neurosurgery Wiki**

Permanent link:

[https://neurosurgerywiki.com/wiki/doku.php?id=pulsed\\_arterial\\_spin\\_labeling](https://neurosurgerywiki.com/wiki/doku.php?id=pulsed_arterial_spin_labeling)

Last update: **2024/06/07 02:52**

