

Silent Scan (GE Healthcare, Milwaukee, Wisconsin) [MRA \(Silent MRA\)](#) is a non-contrast-enhanced MRA technique that uses an [ultrashort TE](#) (UTE) combined with [arterial spin labeling](#), which is used as a preparation pulse to visualize [blood flow](#). Subsequent data acquisition is based on 3D radial sampling.

From:

<https://neurosurgerywiki.com/wiki/> - **Neurosurgery Wiki**

Permanent link:

https://neurosurgerywiki.com/wiki/doku.php?id=ultrashort_te

Last update: **2024/06/07 02:57**

