

The next generation in dural repair is now FDA cleared! Introducing Cerafix® Dura Substitute The FIRST synthetic, resorbable, electrospun dura substitute.

Synthetic: 100% resorbable synthetic material engineered to mimic native dura<sup>2</sup>

Soft: Soft, compliant, drapable material offers procedural versatility; conforms to complex defects<sup>2</sup>

Strong: Sufficient tensile strength and tear resistance for suture application<sup>2</sup>

Safe: No risk of allergic reaction or disease transmission associated with biologics<sup>2 3</sup>

Facilitates cellular infiltration, tissue integration, and neoduralization.<sup>2 3</sup>

Product Code Product Size

C1-1×1 1in x 1in (2.5cm x 2.5cm)

C1-1×3 1in x 3in (2.5cm x 7.5cm)

C1-2×2 2in x 2in (5.0cm x 5.0cm)

C1-3×3 3in x 3in (7.5cm x 7.5cm)

C1-4×5 4in x 5in (10.0cm x 12.5cm)

From:

<https://neurosurgerywiki.com/wiki/> - **Neurosurgery Wiki**

Permanent link:

<https://neurosurgerywiki.com/wiki/doku.php?id=cerafix>

Last update: **2024/06/07 02:53**

