

Maximum Cranial Length (g-op): distance between glabella (g) and opisthocranium (op) in the midsagittal plane, measured in a straight line. Instrument: spreading caliper. Comment: Place skull on side, holding one end of caliper at glabella and extending caliper until maximum diameter at posterior aspect of skull is obtained.

From:

<https://neurosurgerywiki.com/wiki/> - **Neurosurgery Wiki**

Permanent link:

https://neurosurgerywiki.com/wiki/doku.php?id=cranial_length

Last update: **2024/06/07 02:51**

