



For three-dimensional fixation of [spatulas](#)

The position of the flexible holding arm and the LEYLA retractor can be adjusted in three dimensions. A locking mechanism allows fixating the devices in any position.

For endoscopes and MINOP® system, ventriculoscapy or endoscope-assisted neurosurgery

Maximum ease of handling

Adjustability of clamping and holding force

For spatula fixation

Various clamping elements available

From:

<https://neurosurgerywiki.com/wiki/> - **Neurosurgery Wiki**

Permanent link:

https://neurosurgerywiki.com/wiki/doku.php?id=leyla_retractor

Last update: **2024/06/07 02:48**

